

HEMATOLOGY ANALYZER ANA11-35T



HEMATOLOGY ANALYZER ANA11-35T

Blood Analyzer is used to determine components in blood. It is used by hospitals, medical labs, forensic labs and used to test for many things, such as blood cell counts, therapeutic drug monitoring, illegal drug use, blood typing, protein analysis, checking thyroid function, checking for the presence of antibodies.

Used in Forensic Lab, Protein Analysis, Hospitals, Medical, Laboratory, Research, Industry, Pharmaceutical.

Also known as Laboratory Blood Analyzer, Laboratory Hematology Analyzer.

ANA11-35T HEMATOLOGY ANALYZER

High Efficiency and convenience: (1) 60 samples per hour, minimize your workload. (2) 1 click for analysis with 20 parameters+3 histograms.

Complete diagnostic solution: (1) Two reagents only(500mL LYSE, 20 L DILUENT), cost-effective. (2) International brand control and calibrator provided, high quality.

Powerful data management: (1) Store 100,000 results, large memory. (2) Easy to check history results, humanized design. (3) Support LIS, easy data transmission.

Intelligent design: (1) Intelligent software menu for hardware diagnosis. (2) Advanced high voltage and flush to remove blockage automatically.



SPECIFICATIONS

Model	ANA11-35T
Interface	USB port, RS232
Test Principle	Impedance for cell counting and cyanide-free method for HGB
Clog Clear	High voltage and high pressure flush
Parameter	WBC, Neu#, Lym#, Mid#, Neu%,Lym%, Mid%, RBC, HGB, HCT, MCV, MCH, MCHC,RDW-SD, RDW-CV, PLT, MPV, PDW, PCT, P-LCR
Histograms	WBC, RBC, PLT
Sample mode	Venous mode: 10pL, venous blood Capillary mode: 10pL, capillary blood Prediluted mode: 20pL capillary blood
Language	English, Spanish, Portuguese, Italian, etc.
QC Mode	3 level QC, L-J graph
Calibration Mode	Auto mode and Manual mode
Reagent	Lyse(500ml) and Diluent(20L) Reagent required for single test: Diluent: 26mL, Lyse: 0.35mL
Memory	100,000 sample results
Display	Touch screen LCD
Printer	Thermal printer, supports external printer
Dimension	514x594x645 mm
Weight	24.7 kg
Power Supply	100~240 AC, 50/60±1 Hz

ACCESSORIES

Accessory Code	Name	Unit	1
1500705017	Power Cable (1.8 m)	pc	1

1500705018	Ground Cable	pc	1
1500705019	Serial Port Cable	pc	1
1500705020	Mouse	pc	1
1500705021	Keyboard	pc	1
1500705022	Touch Pen	pc	1
1500705023	Side Door Key	pc	1
1500705024	Hex Wrench (Φ 2.5)	set	1
1500705025	Service Pack	set	1
1500705026	Thermal Printer Paper (Width: 80 mm)	roll	1
1500705027	Diluent (20 L)	bottle	1
1500705028	Lyse Solution (500 mL)	bottle	1
1500705029	PB Cleanser (50 mL)	bottle	1
1500705030	E-Z Cleanser (50 mL)	bottle	1
1500705031	Diluent Cap Component	set	1
1500705032	Lyse Solution Cap Component	set	1
1500705033	Waste Bottle Cap Component	set	1
1500705034	Waste Bottle	pc	1

OPTIONAL ACCESSORIES

Accessory Code	Name	Unit
1500705005	Thermal Printer Paper (Width: 80 mm)	roll
1500705006	PB Cleanser (50 mL)	bottle
1500705007	E-Z Cleanser (50 mL)	bottle
1500705008	Service Manual	pc
1500705009	Diluent (20 L)	bottle
1500705010	Lyse solution (500 mL)	bottle

Accessory Code	Name	Unit
1500705011	CBC-3D Blood QC High Level (3.0 mL) (R&D)	tube
1500705012	CBC-3D Blood QC Normal Level (3.0 mL) (R&D)	tube
1500705013	CBC-3D Blood QC Low Level (3.0 mL) (R&D)	tube
1500705014	Lyse solution (1 L)	bottle
1500705015	CBC-CAL PLUS Calibrator (3.0 mL) (R&D)	tube
1500705016	Barcode Scanner Assembly	pc

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@m.labstac.com | Website: labstac.com