

# VERTICAL AUTOCLAVE AUT033-050



# VERTICAL AUTOCLAVE AUT033-050

## AUT033-050 VERTICAL AUTOCLAVE

Three times pre-vacuum  
High quality stainless steel material  
Microcomputer automatic control, arbitrary setting sterilization parameters  
Safety interlocking device  
LED touch screen display  
Standard test interface  
Self-expanding seal  
With drying function  
With stainless steel bucket or basket  
With automatic protection function over temperature protection; over pressure self discharge protection ; low water level protection; anti dry burning.  
Automatic water  
External steam Generator  
Buzzer reminder after sterilization, automatic stop.  
Automatic discharge of cold air, automatic steam exhaust after sterilization  
Internal printer



## SPECIFICATIONS

Model	AUT033-050
Capacity	50 L
Power	3.5 kW
Voltage	AC 220V. 50HZ
Design pressure	0.25 Mpa
Design temperature	139°C
Rated working pressure	0.22 Mpa
Rated working temperature	134°C
Sterilization temperature setting range	116°C-134°C
Sterilization time setting range	4 ~ 120 min
Drying time setting range	0 ~ 240 min
Chamber volume	Ø 400x505
Bucket dimension	Ø 380x430
Basket dimension	Ø 360x280x1
Outer dimension	630x500x1180
Packing size	800x600x1230
Net weight	118 Kg
Gross Weight	138 Kg

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)