

HORIZONTAL AUTOCLAVE AUT22-450



HORIZONTAL AUTOCLAVE AUT22-450

Laboratory horizontal autoclaves are useful in bulk sterilization applications of equipment and materials in various industries and research institutions. Also for disposing biological waste, reducing the risk of contamination and sterile the waste before discarding. We provide a variety of horizontal autoclave according to applications, features and customers requirements.

Used in Bulk sterilization, Microbiology Lab, Research Institutions, Laboratory, Pharmaceutical, Food industry, Chemical industry.

Also known as Laboratory Sterilizer, Horizontal Pressure Sterilizer, Laboratory Horizontal Autoclave,.

AUT22-450 HORIZONTAL AUTOCLAVE



SPECIFICATIONS

Model	AUT22-450
Capacity	450 L
Chamber Size (Wx H x L)	672x672x1000 mm
External Size (WxHxL)	1314x1260x1980 mm
Max Designed Pressure	Inner chamber - 0.1/0.3 Mpa, Jacket 0.3MPa
Working Pressure	0.23 Mpa
Built - in Steam Generator	0.8 Mpa
Designed Temperature	150 °C
Maximum Working Temperature	134 °C
Temperature Control Precision	± 0.2 °C
Consumption	35kVA
Standard Accessories	Trolley, Printer
Altitude	≤ 2000m; could be customised @ 2000-3000m
Power Supply	380 V 50 Hz
Package Size	Main Instrument: 1424x1460x2295; Trolley: 1360x640x1245
Gross Weight	Main Instrument: 1220; Trolley: 110

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com