

CLASS B MEDICAL AUTOCLAVE AUT32-08B



CLASS B MEDICAL AUTOCLAVE AUT32-08B

Medical Autoclave assure proper cleaning and disinfection of instruments and utensils for medical, dental and laboratory field. It minimizes the risk of exposure to infection for each patient and provider. Due to its door locking system and small utilization of space, it becomes an ideal choice for the customers.

Used in Medical, Dental, Laboratory, Pharmaceutical, Research.

Also known as Medical Sterilizer, Dental Autoclave, Medical Autoclave, Pre and Post Vacuum Autoclave, Benchtop Autoclave, Benchtop Sterilizer, Benchtop Steam Sterilizer, Laboratory Medical Autoclave.

AUT32-08B CLASS B MEDICAL AUTOCLAVE

Pre and Post Vacuum Sterilizer

Dual core control system

11 Preset Sterilization Cycles

Maximum negative pressure of -0.80 bar

Preheated chamber keeps the chamber warm and ready to use

LCD display

High density bacteria filter

User-defined cycles and cleaning programs

Drying time: 1-60 min

Printer: Thermal/Stylus

Interface: RS232

ISO & CE 1366 Certified

Vacuum Helix and Bowie & Dick Test passed

Powerful steam generator and vacuum pump

First class safety system with Electrical lock

Multi languages available: English, Italian, German, Chinese

Separate and set up and service menu

Password protection

Safety features: Over heat control, over pressure control, safety door locking system, low water level indicator



SPECIFICATIONS

Model	AUT32-08B
Capacity	8 L
Sterilization Temperature	121°C, 134°C
Sterilization Time	4 min - 60 min
Working Pressure	0.21 - 0.23 MPa
Water Tank	2 x 4.5L Internal Stainless Steel water tanks
USB output	yes
Lock system of door	Mechanic + Electrical lockers
Chamber Dimension	170x320 mm
Overall Dimension (LxWxH)	535x365x320 mm
Packing Size	780 x 500 x 430 mm
Chamber Material	Stainless Steel SUS 304
Weight	59/70 kg

Power	2400 W
Power Supply	220 V, 50Hz

OPTIONAL ACCESSORIES

Accessory Code	Name	Unit
1600505005	Ring	1
1600505006	Inner temperature sensor	1
1600505007	pt1000	1
1600505008	pt1000 (packing with plastic)	1
1600505009	Heating circle	1
1600505010	Steam generator cover	1Set
1600505011	Steam generator stand	1Set
1600505012	Steam generator spring type gasket	1
1600505013	Heating Rod	3
1600505014	Water Pump	1
1600505015	Vacuum Pump	1
1600505016	Fuse Wire	2
1600505017	Normally closed solenoid valve	3

Accessory Code	Name	Unit
1600505018	Three-way solenoid valve	2
1600505019	Germ-tight filter	1
1600505020	Voltage transformer	1
1600505021	switch	1
1600505022	Bleeder valve	2
1600505023	Fan	1
1600505024	Printer	1
1600505025	Control panel	1set
1600505026	LED	1set
1600505027	Big Pipe	1Kg
1600505028	Small Pipe	1kg
1600505029	Liquid sensor	3
1600505030	Electromagnet	1

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com