

# CLASS N MEDICAL AUTOCLAVE AUT35-003



# CLASS N MEDICAL AUTOCLAVE AUT35-003

Medical Autoclave assure proper cleaning and disinfection of instruments and utensils for medical, dental and laboratory field. It minimizes the risk of exposure to infection for each patient and provider. Due to its door locking system and small utilization of space, it becomes an ideal choice for the customers.

Used in Medical, Dental, Laboratory, Pharmaceutical, Research.

Also known as Medical Sterilizer, Dental Autoclave, Medical Autoclave, Pre and Post Vacuum Autoclave, Benchtop Autoclave, Benchtop Sterilizer, Benchtop Steam Sterilizer, Laboratory Medical Autoclave.

## AUT35-003 CLASS N MEDICAL AUTOCLAVE



Product Image Coming Soon

### SPECIFICATIONS

|                           |                                   |
|---------------------------|-----------------------------------|
| Model                     | AUT35-003                         |
| Capacity                  | 3 L                               |
| Sterilization Pressure    | 0.22 Mp                           |
| Sterilization Temperature | 134 °C                            |
| Sterilization Time        | 4 min                             |
| Sterilization Cycle time  | 11 min                            |
| Required water            | Distilled water or purified water |
| Product Size              | 300X190X350 mm                    |
| Outer Dimension           | 370x260x460 mm                    |
| Gross Weight              | 7.5 Kg                            |
| Power                     | 0.5 KW                            |
| Voltage                   | AC220/110V+ 10%                   |

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)