

CLASS N MEDICAL AUTOCLAVE AUT35-23N



CLASS N MEDICAL AUTOCLAVE AUT35-23N

Medical Autoclave assure proper cleaning and disinfection of instruments and utensils for medical, dental and laboratory field. It minimizes the risk of exposure to infection for each patient and provider. Due to its door locking system and small utilization of space, it becomes an ideal choice for the customers.

Used in Medical, Dental, Laboratory, Pharmaceutical, Research.

Also known as Medical Sterilizer, Dental Autoclave, Medical Autoclave, Pre and Post Vacuum Autoclave, Benchtop Autoclave, Benchtop Sterilizer, Benchtop Steam Sterilizer, Laboratory Medical Autoclave.

AUT35-23N CLASS N MEDICAL AUTOCLAVE



Product Image Coming Soon

SPECIFICATIONS

Model	AUT35-23N
Capacity	23 L
Product Size	450x550x350 mm
Chamber Dimension	φ249x450 mm
Package Size	670x550x500 mm
Gross Weight	33 Kg
Net Weight	26 Kg
Power	800 W
Voltage	AC220/110V+10%, 50/60 Hz+ 1Hz

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com