

PORTABLE AUTOCLAVE AUT43-18



PORTABLE AUTOCLAVE AUT43-18

Portable Autoclave is used for general purpose sterilization. Saturated steam is used to sterilize articles rapidly and efficiently. It is also used as a cooker in plateau areas and also for making high quality drinking water in Drinking Water Plant.

Used in Medical, Research Institutions, Laboratory, Drinking Water Plant, Food industry, Hospital, Agriculture, Hospital, Scientific Research..

Also known as Portable Sterilizer, Steam Autoclave Sterilizer, Portable Autoclave, Pressure Steam Sterilizer, Portable Pressure Steam Sterilizer, Dental Autoclave.

AUT43-18 PORTABLE AUTOCLAVE

SPECIFICATIONS

| | |
|------------------------------|-----------------|
| Model | AUT43-18 |
| Capacity | 18 L |
| Function | Electric heated |
| Sterilization Temperature | 126°C |
| Steam Exhausting Temperature | ≤1°C |
| Working Pressure | 0.14-0.16Mpa |
| Maximum Working Pressure | 0.165 Mpa |
| Timer Range | 0-99min |
| Chamber Dimension | 280x260mm |
| Overall Dimension (LxWxH) | 440x420x630 mm |
| Chamber Material | Stainless Steel |
| Weight | 12/14 kg |
| Power | 2 KW |
| Power Supply | 220 V, 50 Hz |

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com