

# LOW TEMPERATURE CIRCULATOR BAC027-23



# LOW TEMPERATURE CIRCULATOR BAC027-23

## BAC027-23 LOW TEMPERATURE CIRCULATOR

Microprocessor temperature controller with LCD screen ensures precise and reliable control, easy to operate

No angle in bath chamber, easy to clean

With interface to external water bath

Easy to set adjustable timer. (1 minute to 5,999 minutes)

Maintenance-free operation with easy to clean surface

Both heater and bath chamber are made of corrosion-resistant stainless steel

Drain valve makes emptying of bath water fast and easy for cleaning and moving

Programming setting function with 7 periods and 9 steps for each period, which means there are 63 programmable steps in total



## SPECIFICATIONS

Model	BAC027-23
Capacity	22.5 L
Temperature range	RT+5 - 150°C
Precision	±0.1
Liquid tank opening depth	330x300/150 mm
Power requirements	AC 220 V 50 HZ
Pump	8 L/min
Power Consumption	1050 W

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)