

WATER BATH BAC028-12



WATER BATH BAC028-12

BAC028-12 WATER BATH

PID temperature controller, large LCD display screen and interface provides for userfriendly operation

Easy to set Adjustable Timer. (1 minute to 5,999 minutes)

Maintenance-free operation with easy to clean surface

Both heater and bath chamber are made of corrosion-resistant stainless steel

No angle in bath chamber, easy to clean

The standard stainless-steel bottom plate, helps prevents direct contact by accessories and tubes to heating element and sensors

Drain valve makes emptying of bath water fast and easy for cleaning and moving



SPECIFICATIONS

Model	BAC028-12
Capacity	12 L
Power requirements	AC 220 V 50 HZ
Power	800 W
Temperature Range	RT+5 - 100°C
Temperature fluctuation	±0.2°C
Temperature resolution	0.1°C
Internal dimension	305Wx150Dx240H mm
External dimension	353Wx340Dx265H mm
Timer	1 - 5999 min
Net weight	12 Kg
Porous cover(lid)	4 - hole

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com