

LOW TEMPERATURE CIRCULATOR BAC032-05A



LOW TEMPERATURE CIRCULATOR BAC032-05A

BAC032-05A LOW TEMPERATURE CIRCULATOR

Microprocessor temperature controller

Audible and visible alarm for temperature and water level

R134a refrigerant, imported compressor

With interface to external water bath

RS485 connector is option which can connect computer to record the parameters and the variations of temperature



SPECIFICATIONS

Model	BAC032-05A
Capacity	4.5 L
Temperature range	-10-100°C
Precision	±0.2
Interior Dimension	150x160x150
Electrical Requirements	220 V 50 Hz
Pump (flux)	8 L/min
Power Consumption	2300 W

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@mail.labstac.com | Website: labstac.com