

# LOW TEMPERATURE CIRCULATOR BAC032-05D



# LOW TEMPERATURE CIRCULATOR BAC032-05D

## BAC032-05D LOW TEMPERATURE CIRCULATOR

Microprocessor temperature controller

Audible and visible alarm for temperature and water level

R134a refrigerant, imported compressor

With interface to external water bath

RS485 connector is option which can connect computer to record the parameters and the variations of temperature



## SPECIFICATIONS

Model	BAC032-05D
Capacity	4.5 L
Temperature range	-40-100°C
Precision	±0.2
Interior Dimension	150x160x150
Electrical Requirements	220 V 50 Hz
Pump (flux)	8 L/min
Power Consumption	3150 W

# LABSTAC

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)