

LOW TEMPERATURE CIRCULATOR BAC032-13D



LOW TEMPERATURE CIRCULATOR BAC032-13D

BAC032-13D LOW TEMPERATURE CIRCULATOR

Microprocessor temperature controller

Audible and visible alarm for temperature and water level

R134a refrigerant, imported compressor

With interface to external water bath

RS485 connector is option which can connect computer to record the parameters and the variations of temperature



SPECIFICATIONS

Model	BAC032-13D
Capacity	13 L
Temperature range	-50-100°C
Precision	±0.2
Interior Dimension	240x170x200
Electrical Requirements	220 V 50 Hz
Pump (flux)	8 L/min
Power Consumption	3100 W

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@mail.labstac.com | Website: labstac.com