

LOW TEMPERATURE CIRCULATOR BAC032-15B



LOW TEMPERATURE CIRCULATOR BAC032-15B

BAC032-15B LOW TEMPERATURE CIRCULATOR

Microprocessor temperature controller

Audible and visible alarm for temperature and water level

R134a refrigerant, imported compressor

With interface to external water bath

RS485 connector is option which can connect computer to record the parameters and the variations of temperature



SPECIFICATIONS

Model	BAC032-15B
Capacity	14.5 L
Temperature range	RT+5-100°C
Precision	±0.2
Interior Dimension	240x170x200
Electrical Requirements	220 V 50 Hz
Pump (flux)	8 L/min
Power Consumption	1050 W

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@mail.labstac.com | Website: labstac.com