

# LOW TEMPERATURE CIRCULATOR BAC033-13B



# LOW TEMPERATURE CIRCULATOR BAC033-13B

## BAC033-13B LOW TEMPERATURE CIRCULATOR

High precision temperature controller



### SPECIFICATIONS

Model	BAC033-13B
Capacity	13 L
Temperature range	-30-100°C
Precision	±0.02
Peristome Dimension	240x170x200
Pump(flux)	15 L/min
Power	2850 W

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)