

LOW TEMPERATURE CIRCULATOR BAC033-13C



LOW TEMPERATURE CIRCULATOR BAC033-13C

BAC033-13C LOW TEMPERATURE CIRCULATOR

High precision temperature controller



SPECIFICATIONS

Model	BAC033-13C
Capacity	13 L
Temperature range	-40-100°C
Precision	±0.02
Peristome Dimension	240x170x200
Pump(flux)	15 L/min
Power	2850 W

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com