

SHAKING WATER BATH BAC034-33A



SHAKING WATER BATH BAC034-33A

BAC034-33A SHAKING WATER BATH

A stamping molding stainless steel tank, easy to clean

LCD screen, multiple data display with timing function, easy to operate

Stainless steel shelves cover heater and sensor to avoid damage during using

Once-forming stainless steel lid

Cut off heater automatically in case of lack of water, meanwhile visible and audible alarm ensures to remind users in time

Independent temperature-limiting alarm system

temperature error alarm

Test tube holder can be placed



SPECIFICATIONS

Model	BAC034-33A
Capacity	33 L
Temperature range	RT+5-100°C
Display Resolution	0.1°C
Temperature Uniformity	±1°C
Shaking speed range	30-150rpm
Amplitude	30mm (Standard) or 40mm (Option)
Power Consumption	1250 W
Interior Dimension	438Wx310Dx250H mm
Exterior Dimension	643Wx350Dx353H mm

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com