

SHAKING WATER BATH BAC034-47A



SHAKING WATER BATH BAC034-47A

BAC034-47A SHAKING WATER BATH

- A stamping molding stainless steel tank, easy to clean
- LCD screen, multiple data display with timing function, easy to operate
- Stainless steel shelves cover heater and sensor to avoid damage during using
- Once-forming stainless steel lid
- Cut off heater automatically in case of lack of water, meanwhile visible and audible alarm ensures to remind users in time
- Independent temperature-limiting alarm system
- temperature error alarm
- Test tube holder can be placed



SPECIFICATIONS

Model	BAC034-47A
Capacity	47 L
Temperature range	RT+5-100°C
Display Resolution	0.1°C
Temperature Uniformity	±1°C
Shaking speed range	30-150rpm
Amplitude	30mm (Standard) or 40mm (Option)
Power Consumption	1650 W
Interior Dimension	618Wx310Dx250H mm
Exterior Dimension	823Wx350Dx355H mm

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com