

SHAKING WATER BATH BAC034-47B



SHAKING WATER BATH BAC034-47B

BAC034-47B SHAKING WATER BATH

A stamping molding stainless steel tank, easy to clean

LCD screen, multiple data display with timing function, easy to operate

Stainless steel shelves cover heater and sensor to avoid damage during using

Once-forming stainless steel lid

Cut off heater automatically in case of lack of water, meanwhile visible and audible alarm ensures to remind users in time

Independent temperature-limiting alarm system

temperature error alarm

Test tube holder can be placed



SPECIFICATIONS

Model	BAC034-47B
Capacity	47 L
Temperature range	RT+5-100°C
Display Resolution	0.1°C
Temperature Uniformity	±1°C
Shaking speed range	30-150rpm
Amplitude	30mm (Standard) or 40mm (Option)
Power Consumption	1650 W
Interior Dimension	618Wx310Dx250H mm
Exterior Dimension	823Wx350Dx355H mm

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@mail.labstac.com | Website: labstac.com