

WATER BATH BAC035-05A



WATER BATH BAC035-05A

BAC035-05A WATER BATH

- A stamping molding stainless steel tank, easy to clean
- LCD screen, multiple data display with timing function, easy to operate
- Stainless steel shelves cover heater and sensor to avoid damage during using
- Once-forming stainless steel lid
- Cut off heater automatically in case of lack of water, meanwhile visible and audible alarm ensures to remind users in time
- Independent temperature-limiting alarm system
- temperature error alarm
- Test tube holder can be placed

SPECIFICATIONS

Model	BAC035-05A
Capacity	2 holes Φ 112 mm
Temperature range	RT+5-100°C
Electrical Requirement	AC 220 50 Hz
Power Consumption	500 W
Temperature Stability	\pm 0.3°C
Temp Alarm	\pm 2°C
Interior dimension	130Wx280Dx150H mm
Exterior dimension	396Wx250Dx260H mm
Timing Range	1~5999 min

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com