

WATER BATH BAC035-30



WATER BATH BAC035-30

BAC035-30 WATER BATH

A stamping molding stainless steel tank, easy to clean

LCD screen, multiple data display with timing function, easy to operate

Stainless steel shelves cover heater and sensor to avoid damage during using

Once-forming stainless steel lid

Cut off heater automatically in case of lack of water, meanwhile visible and audible alarm ensures to remind users in time

Independent temperature-limiting alarm system

temperature error alarm

Test tube holder can be placed



SPECIFICATIONS

Model	BAC035-30
Capacity	20 L
Temperature range	RT+5-100°C
Electrical Requirement	AC 220 50 Hz
Power Consumption	1000 W
Temperature Stability	±0.2°C
Temp Alarm	0.1°C
Interior dimension	500Wx300Dx200H mm
Exterior dimension	680Wx360Dx390H mm
Timing Range	1~5999 min

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com