

WATER BATH BAC035-30



WATER BATH BAC035-30

BAC035-30 WATER BATH

A stamping molding stainless steel tank, easy to clean

LCD screen, multiple data display with timing function, easy to operate

Stainless steel shelves cover heater and sensor to avoid damage during using

Once-forming stainless steel lid

Cut off heater automatically in case of lack of water, meanwhile visible and audible alarm ensures to remind users in time

Independent temperature-limiting alarm system

temperature error alarm

Test tube holder can be placed



SPECIFICATIONS

| | |
|------------------------|-------------------|
| Model | BAC035-30 |
| Capacity | 20 L |
| Temperature range | RT+5-100°C |
| Electrical Requirement | AC 220 50 Hz |
| Power Consumption | 1000 W |
| Temperature Stability | ±0.2°C |
| Temp Alarm | 0.1°C |
| Interior dimension | 500Wx300Dx200H mm |
| Exterior dimension | 680Wx360Dx390H mm |
| Timing Range | 1~5999 min |

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com