

WATER BATH BAC035-69



WATER BATH BAC035-69

BAC035-69 WATER BATH

- A stamping molding stainless steel tank, easy to clean
- LCD screen, multiple data display with timing function, easy to operate
- Stainless steel shelves cover heater and sensor to avoid damage during using
- Once-forming stainless steel lid
- Cut off heater automatically in case of lack of water, meanwhile visible and audible alarm ensures to remind users in time
- Independent temperature-limiting alarm system
- temperature error alarm
- Test tube holder can be placed



SPECIFICATIONS

Model	BAC035-69
Capacity	11 L
Temperature range	RT+5-80°C
Electrical Requirement	AC 220 50 Hz
Power Consumption	800 W
Temperature Stability	±0.2°C
Temp Alarm	0.1°C
Interior dimension	300Wx240Dx200H mm
Exterior dimension	480Wx300Dx480H mm
Timing Range	1~5999 min

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com