

WATER BATH BAC036-06



WATER BATH BAC036-06

BAC036-06 WATER BATH

- Microprocessor control with timing function
- Digital display
- Audible and visible alarm for over temperature



SPECIFICATIONS

Model	BAC036-06
Capacity	2.1 L x 3
Electrical Requirement	220 V 50 HZ
Power Consumption	750 W
Temperature Range	RT+5-99°C
Temperature Stability	±0.3°C
Display Resolution	0.1 °C
Interior dimension	150Wx125Dx110H mm
Exterior dimension	490Wx245Dx310H mm
Timing Range	1~999 min

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com