

# WATER BATH BAC036-11



# WATER BATH BAC036-11

## BAC036-11 WATER BATH

Microprocessor control with timing function  
Digital display  
Audible and visible alarm for over temperature



### SPECIFICATIONS

Model	BAC036-11
Capacity	11 L
Electrical Requirement	220 V 50 HZ
Power Consumption	500 W
Temperature Range	RT+5-100°C
Temperature Stability	±0.5°C
Display Resolution	0.1 °C
Interior dimension	420Wx180Dx150H mm
Exterior dimension	570Wx220Dx260H mm
Timing Range	1~999 min

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)