

# WATER BATH BAC036-22



# WATER BATH BAC036-22

## BAC036-22 WATER BATH

- Microprocessor control with timing function
- Digital display
- Audible and visible alarm for over temperature



### SPECIFICATIONS

Model	BAC036-22
Capacity	22 L
Electrical Requirement	220 V 50 HZ
Power Consumption	600 W
Temperature Range	RT+5-100°C
Temperature Stability	±0.5°C
Display Resolution	0.1 °C
Interior dimension	450Wx300Dx190H mm
Exterior dimension	610Wx340Dx260H mm
Timing Range	1~999 min



**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)