

WATER BATH BAC036-30



WATER BATH BAC036-30

BAC036-30 WATER BATH

- Microprocessor control with timing function
- Digital display
- Audible and visible alarm for over temperature



SPECIFICATIONS

Model	BAC036-30
Capacity	30 L
Electrical Requirement	220 V 50 HZ
Power Consumption	1000 W
Temperature Range	RT+5-100°C
Temperature Stability	±0.5°C
Display Resolution	0.1 °C
Interior dimension	600Wx300Dx190H mm
Exterior dimension	760Wx340Dx280H mm
Timing Range	1~999 min



Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com