

# WATER BATH BAC036-34



# WATER BATH BAC036-34

## BAC036-34 WATER BATH

- Microprocessor control with timing function
- Digital display
- Audible and visible alarm for over temperature



### SPECIFICATIONS

|                        |                   |
|------------------------|-------------------|
| Model                  | BAC036-34         |
| Capacity               | 34 L              |
| Electrical Requirement | 220 V 50 HZ       |
| Power Consumption      | 1000 W            |
| Temperature Range      | RT+5-100°C        |
| Temperature Stability  | ±0.5°C            |
| Display Resolution     | 0.1 °C            |
| Interior dimension     | 600Wx300Dx190H mm |
| Exterior dimension     | 750Wx350Dx300H mm |
| Timing Range           | 1~999 min         |



**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)