

HIGH TEMPERATURE CIRCULATOR BAC037-22



HIGH TEMPERATURE CIRCULATOR BAC037-22

BAC037-22 HIGH TEMPERATURE CIRCULATOR

- Microprocessor control with timing function
- Digital display
- Audible and visible alarm for over temperature



SPECIFICATIONS

Model	BAC037-22
Capacity	22 L
Electrical Requirement	220 V 50 HZ
Power Consumption	1000 W
Temperature Range	RT+5-70°C
Temperature Stability	±0.3°C
Display Resolution	0.1 °C
Interior dimension	500Wx300Dx150H mm
Exterior dimension	660Wx340Dx200H mm
Timing Range	1~999 min

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com