

# HIGH TEMPERATURE CIRCULATOR BAC037-30



# HIGH TEMPERATURE CIRCULATOR BAC037-30

## BAC037-30 HIGH TEMPERATURE CIRCULATOR

Microprocessor control with timing function  
Digital display  
Audible and visible alarm for over temperature



### SPECIFICATIONS

Model	BAC037-30
Capacity	30 L
Electrical Requirement	220 V 50 HZ
Power Consumption	1000 W
Temperature Range	RT+5-70°C
Temperature Stability	±0.3°C
Display Resolution	0.1 °C
Interior dimension	600Wx300Dx190H mm
Exterior dimension	770Wx370Dx280H mm
Timing Range	1~999 min

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)