

# DRY BATH CHILLER BAC71-00



# DRY BATH CHILLER BAC71-00

Dry Bath Chillers are reliable, easy to use and maintain. Specially designed to meet application demands of laboratories. Greater energy efficiency and less environmental impact, as well as offering increased temperature ranges and greater stability.

Used in Pharmaceutical, Chemical, Food Safety, Environment, Quality inspection, Laboratory.

Also known as Laboratory Dual Temperature Dry Bath Incubator.

## BAC71-00 DRY BATH CHILLER



### SPECIFICATIONS

Model	BAC71-00
Capacity	4 L
Function	Basic model
Temperature Range	Maintains at 0°C
Temperature Uniformity	±0.5°C
Cleaning Tank (LxWxH)	300x150x100 mm
Outer Dimension (WxDxH)	420x230x270 mm
Weight	1.4 kg
Temperature Precision	±0.5°C
Operation Ambient Temperature	8-30°C
Maximum Input Power	220
Maximum Output Power	180
Heating Power	NO
Power Adjustable (%)	NO
Drain	NO
Level Control System	NO
Power Supply	12 V DC, 7A



**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)