

ANALYTICAL BALANCE BAL12-220IC



ANALYTICAL BALANCE BAL12-220IC

Analytical Balances are designed to meet the specific weighing requirement in the laboratory working environment. Performance, durability, robust, easy-to-use, precise results, efficient handling, space saving makes analytical balance an ideal choice.

Used in Totalization, Formulation, SQC, Industry, Laboratory, Research.

Also known as Lab Balance, Laboratory Analytical Balance, Analytical Scale, Laboratory Analytical Scale.

BAL12-220IC ANALYTICAL BALANCE

Electromagnetic balanced sensors

The multiple functions modes with tare, counting, unit conversion and weighing

External Calibration



SPECIFICATIONS

Model	BAL12-220IC
Capacity (g)	220 g
Resolution (g)	0.0001 g
Min Weighing (g)	0.0004 g
Stable Time	≤3 sec
Display	LCD(white back light with black font)
Pan Size	Φ 80 mm
Operation Temperature	5-25°C
Repeatability	±0.0002 g
Liner	±0.0003 g
Cal. Weight	Internal Calibration
Overall Dimension	240x190x265 mm
Packing Size	Inside: 415x285x430 mm Outside: 475x315x485 mm
Gross Weight	10 kg

OPTIONAL ACCESSORIES

Accessory Code	Name
3500505005	RS232C & RS485 Interface
3500505006	Printer

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com