

**LABSTAC**

# PRECISION BALANCE BAL24-1100



labstac.com

# PRECISION BALANCE BAL24-1100

Precision Balance incorporates unique high end technology which is built to match the highest international standards. Outstanding weighing performance and innovative sample handling technologies are used. Due to its efficiency, error-free, and compliant features it becomes an unique product.

Used in Chemistry, Biological Research, Industry, Laboratory.

Also known as Loading Balance, Laboratory Precision Balance, Laboratory Loading Balance, Precision Scale, Laboratory Precision Scale.

## BAL24-1100 PRECISION BALANCE

Stainless steel Pan

LCD digital display

The multiple functions modes with tare, counting, unit conversion and weighing

4 conversion units: oz,ct,lb,g



## SPECIFICATIONS

Model	BAL24-1100
Capacity (g)	1100 g
Resolution (g)	0.01 g
Min Weighing (g)	0.04 g
Stable Time	≤2 sec
Pan Size	Φ 130 mm
Operation Temperature	5-25°C
Repeatability	0.02 g
Liner	0.03 g
Cal. Weight	500 g
Packing Size	Inside: 310x230x130 mm Outside: 660x470x320 mm
Gross Weight	2 / 22 kg

## OPTIONAL ACCESSORIES

Accessory Code	Name
3500505005	RS232C & RS485 Interface
3500505006	Printer
3500614005	Windbreak mask

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)