

# PRECISION BALANCE BAL25-0210



# PRECISION BALANCE BAL25-0210

Precision Balance incorporates unique high end technology which is built to match the highest international standards. Outstanding weighing performance and innovative sample handling technologies are used. Due to its efficiency, error-free, and compliant features it becomes an unique product.

Used in Chemistry, Biological Research, Industry, Laboratory.

Also known as Loading Balance, Laboratory Precision Balance, Laboratory Loading Balance, Precision Scale, Laboratory Precision Scale.

## BAL25-0210 PRECISION BALANCE

LCD digital display

The multiple functions modes with tare, counting, unit conversion and weighing



## SPECIFICATIONS

Model	BAL25-0210
Capacity (g)	210 g
Resolution (g)	0.001 g
Min Weighing (g)	0.004 g
Stable Time	≤3 sec
Display	LCD(white back light with black font)
Sensor	Load Cell
Pan Size	Φ 90 mm
Operation Temperature	5-25°C
Repeatability	±0.002 g
Liner	±0.003 g
Cal. Weight	200 g
Overall Dimension	177x173x205 mm
Packing Size	Inside: 350x265x345 mm Outside: 740x550x390 mm
Gross Weight	4 / 16 kg

## OPTIONAL ACCESSORIES

Accessory Code	Name
3500505005	RS232C & RS485 Interface
3500505006	Printer
3500624005	USB

3500624006	Dual Display
------------	--------------



**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)