

# WEIGHING SCALE BAL31-021



# WEIGHING SCALE BAL31-021

Weighing Scale uses high precision weighing sensor to keep weighing stable. Used in laboratory, industrial and food retailing application. Due to its latest design and friendly material, it is widely used in the laboratory. Used in Weighing, Percent Weighing, Industrial, Food Retailing, Laboratory, Research. Also known as Weighing Balance, Laboratory Weighing Scale, Laboratory Weighing Balance.

## BAL31-021 WEIGHING SCALE

Stainless steel Pan

LCD digital display

The multiple functions modes with tare, counting, unit conversion and weighing



## SPECIFICATIONS

Model	BAL31-021
Capacity (g)	21 kg
Resolution (g)	0.1 g
Min Weighing (g)	0.4 g
Stable Time	≤2 sec
Display	LCD(white back light with black font)
Sensor	Load Cell
Pan Size	315x230 mm
Operation Temperature	5-25°C
Repeatability	±0.2 g
Liner	±0.3 g
Overall Dimension	240x190x265 mm
Packing Size	Inside: 410x165x400 mm Outside: 800x420x420 mmm
Gross Weight	4 / 16 kg

## OPTIONAL ACCESSORIES

Accessory Code	Name
3500505005	RS232C & RS485 Interface
3500505006	Printer

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)