

WEIGHING SCALE BAL33-061



WEIGHING SCALE BAL33-061

Weighing Scale uses high precision weighing sensor to keep weighing stable. Used in laboratory, industrial and food retailing application. Due to its latest design and friendly material, it is widely used in the laboratory.

Used in Weighing, Percent Weighing, Industrial, Food Retailing, Laboratory, Research.

Also known as Weighing Balance, Laboratory Weighing Scale, Laboratory Weighing Balance.

BAL33-061 WEIGHING SCALE

Stainless steel Pan

LCD digital display

The multiple functions modes with tare, counting, unit conversion and weighing



SPECIFICATIONS

Model	BAL33-061
Capacity (g)	61 kg
Resolution (g)	10 g
Min Weighing (g)	50 g
Stable Time	≤2 sec
Display	LED
Sensor	Load Cell
Pan Size	410x510 mm / 450x460 mm / 600x800 mm
Operation Temperature	5-25°C
Repeatability	±20 g
Liner	±30 g
Packing Size	420x170x850 mm
Gross Weight	13 kg

OPTIONAL ACCESSORIES

Accessory Code	Name
3500505005	RS232C & RS485 Interface
3500505006	Printer

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com