

WEIGHING SCALE BAL34-160



WEIGHING SCALE BAL34-160

Weighing Scale uses high precision weighing sensor to keep weighing stable. Used in laboratory, industrial and food retailing application. Due to its latest design and friendly material, it is widely used in the laboratory.

Used in Weighing, Percent Weighing, Industrial, Food Retailing, Laboratory, Research.

Also known as Weighing Balance, Laboratory Weighing Scale, Laboratory Weighing Balance.

BAL34-160 WEIGHING SCALE

Stainless steel Pan

LCD digital display

The multiple functions modes with tare, counting, unit conversion and weighing



SPECIFICATIONS

| | |
|-----------------------|----------------|
| Model | BAL34-160 |
| Capacity (g) | 160 kg |
| Resolution (g) | 1 g |
| Min Weighing (g) | 5 g |
| Stable Time | ≤2 sec |
| Display | LED |
| Sensor | Load Cell |
| Pan Size | 410x510 mm |
| Operation Temperature | 5-25°C |
| Repeatability | ±2 g |
| Liner | ±3 g |
| Packing Size | 420x170x850 mm |
| Gross Weight | 13 kg |

OPTIONAL ACCESSORIES

| Accessory Code | Name |
|----------------|--------------------------|
| 3500505005 | RS232C & RS485 Interface |
| 3500505006 | Printer |



Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@labstac.com | Website: labstac.com