

CELL COUNTER BCC31-200



CELL COUNTER BCC31-200

Bright field detection method with 6 chamber slot per slide

Used in Blood analysis, Diagnostic purposes, Research.

Also known as Cell Counter.

BCC31-200 CELL COUNTER

Bright field detection method with 6 chamber slot per slide



SPECIFICATIONS

Model	BCC31-200
Chamber slot	6 slots per slide
Detection method	Brightfield
Focus	Manual focus
Data transfer	USB
Measuring Time per sample	2s
Sample concentration range	1x 10 ⁴ to 3 x 10 ⁷ cells/ml
Particle/Cell diameter range	5-200 µm
Sample volume	25 µ (*6 slots)
Optics	5 MP CMOS
Display	N/A (connected to laptop)
No. of user profiles	More than 100
Weight	1.8 kg

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com