

# CELL COUNTER BCC31-200



# CELL COUNTER BCC31-200

Bright field detection method with 6 chamber slot per slide

Used in Blood analysis, Diagnostic purposes, Research.

Also known as Cell Counter.

## BCC31-200 CELL COUNTER

Bright field detection method with 6 chamber slot per slide

### SPECIFICATIONS

Model	BCC31-200
Chamber slot	6 slots per slide
Detection method	Brightfield
Focus	Manual focus
Data transfer	USB
Measuring Time per sample	2s
Sample concentration range	1x 10 <sup>4</sup> to 3 x 10 <sup>7</sup> cells/ml
Particle/Cell diameter range	5-200 µm
Sample volume	25 µ (*6 slots)
Optics	5 MP CMOS
Display	N/A (connected to laptop)
No. of user profiles	More than 100
Weight	1.8 kg

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)