

LABSTAC

AUTOMATED BLOOD CULTURE SYSTEM BCS11-240



labstac.com

AUTOMATED BLOOD CULTURE SYSTEM

BCS11-240


Developed to accurately detect bacteria, fungus, and other microorganisms in blood and other sterile bodily fluids. For each size laboratory with varying demands, the equipment' numerous capacities may offer dependable, adaptable, and affordable blood culture solutions.

Used in Clinical laboratoy, Medical laboratory.

Also known as Automated Microbial Detection System, Instrumented Blood Culture System, Microbial Detection Instrument Microbial Culture Instrument.

BCS11-240 AUTOMATED BLOOD CULTURE SYSTEM

- Four incubator drawers of 60 cell capacity deliver 240 cell capacity
- One heating system for each incubator drawer to control temperature
- Control module handles up to 64 incubator modules
- Support for random loading and unloading, support for incubation of anonymous bottles, positive anonymous and negative anonymous bottles can be identified
- Double-barcode design to load/unload bottles simply and avoid operation mistakes
- Time to positive (TTP) of more than 90% positive specimens is within 24 hours
- Alarm by audio and light flash



SPECIFICATIONS

Model	BCS11-240
Capacity	240
Dimension (LxWxH)	780x650x980 mm
Weight	200 kg



Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA
Email: contact@mail.labstac.com | Website: labstac.com