

**LABSTAC**

# AUTOMATED BLOOD CULTURE SYSTEM BCS11-60



labstac.com

# AUTOMATED BLOOD CULTURE SYSTEM

## BCS11-60

Developed to accurately detect bacteria, fungus, and other microorganisms in blood and other sterile bodily fluids. For each size laboratory with varying demands, the equipment' numerous capacities may offer dependable, adaptable, and affordable blood culture solutions.

Used in Clinical laboratoy, Medical laboratory.

Also known as Automated Microbial Detection System, Instrumented Blood Culture System, Microbial Detection Instrument Microbial Culture Instrument.

## BCS11-60 AUTOMATED BLOOD CULTURE SYSTEM

60 cell capacity

Continuous detection of the bottles in the cycle of 10 minutes

Swing shaking inoculation is beneficial to the organisms growth

Double-barcode design to load/unload bottles simply and avoid operation mistakes

Automated recognition and alarm of inconsistent bottles

Automated optical detection and manual interpretation to prevent the occurrence of false negative result

Real time dynamic display of the state of culture and growth curve on the graphical user interface

User friendly software system, convenient data management

Rapid selective cultivation of special isolates



### SPECIFICATIONS

Model	BCS11-60
Capacity	60
Dimension (LxWxH)	580x595x590 mm
Weight	65 kg



**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)