

# MICROBIAL ID/AST SYSTEM BMD11-18



# MICROBIAL ID/AST SYSTEM BMD11-18

The device offers quick and precise species-level identification of clinically significant pathogenic bacteria and yeasts based on the most recent microbial categorization system.

Used in Microbial Diagnosis, Epidemiological Investigation, Drug Resistance Trend Analysis and Hospital Infection Management..

## **BMD11-18** MICROBIAL ID/AST SYSTEM

Colorimetric method for microorganism identification and turbidity method for antibiotic susceptibility testing (AST)

Powerful QC function can meet the quality control requirement of any size lab.

Connect with users ♦ LIS and transfer data to WHONET

### **SPECIFICATIONS**

Model	BMD11-18
Detection Method	Photoelectric micro-well by micro-well detection with 8 read detection unit
ID/AST results	Provided in 18-24 hours, as little as 4 to 6 hours

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA  
Email: [contact@mail.labstac.com](mailto:contact@mail.labstac.com) | Website: [labstac.com](http://labstac.com)