

PASS BOX BSC61-600



PASS BOX BSC61-600

It is used to pass goods between room of different level of cleanliness , so as to avoid air convention between them and minimize pollution extent

Used in Cell culturing, Genetics, Microbiological, Research, Cell Biology, Molecular Biology, Plant, Laboratory.

BSC61-600 PASS BOX

Photoelectric sensor, automatic blowing

Adapt with circulation wind to ensure the cleanness of shower area under non-shower status

Double door by electric interlock

LED display with adjustable blowing time(0-99s)

HEPA filter, Class 100 cleanness

SPECIFICATIONS

Model	BSC61-600
UV Lamp	8Wx1
External Size (WxDxH)	760x670x680 mm
Internal Size (WxDxH)	600x600x600 mm
Material	304 stainless steel
Interlock	Electronical
Gross Weight	70kg
Packing Size(WxDxH)	880x800x870 mm
Power Supply	220V, 50Hz; 110V, 50/60Hz

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@m.labstac.com | Website: labstac.com