

# WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET CBT11-12L



# WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET CBT11-12L

Weak Acid and Alkali Chemical Storage Cabinet are specifically designed for the safe storage of Weak Acids, Alkali, Hazardous Chemicals and corrosive materials for example Phosphoric acid, Carbonic acid, Acetic acid, Hydrocyanic acid and other chemical substances.

Used in Laboratory, Research, Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants. Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

## CBT11-12L WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET



### SPECIFICATIONS

Model	CBT11-12L
Capacity	12 gallon (45L)
Material	Double layer anti fire cold-rolled steel PP tray
Adjustable Shelf	1 pc
Door	Single / Manual
External Size	590x460x890 mm
Package Size	710x580x1066 mm
Gross Weight	66 Kg

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)