

# WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET CBT11-30L



# WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET CBT11-30L

Weak Acid and Alkali Chemical Storage Cabinet are specifically designed for the safe storage of Weak Acids, Alkali, Hazardous Chemicals and corrosive materials for example Phosphoric acid, Carbonic acid, Acetic acid, Hydrocyanic acid and other chemical substances.

Used in Laboratory, Research, Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants. Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

## CBT11-30L WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET



### SPECIFICATIONS

Model	CBT11-30L
Capacity	30 gallon (114L)
Material	Double layer anti fire cold-rolled steel PP tray
Adjustable Shelf	1 pc
Door	Double / Manual
External Size	1090x460x1120 mm
Package Size	1210x680x1285 mm
Gross Weight	117 Kg

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)