

WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET CBT11-45L



WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET CBT11-45L

Weak Acid and Alkali Chemical Storage Cabinet are specifically designed for the safe storage of Weak Acids, Alkali, Hazardous Chemicals and corrosive materials for example Phosphoric acid, Carbonic acid, Acetic acid, Hydrocyanic acid and other chemical substances.

Used in Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants, Laboratory, Research. Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

CBT11-45L WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET



SPECIFICATIONS

Model	CBT11-45L
Capacity	45 gallon (170L)
Material	Double layer anti fire cold-rolled steel PP tray
Adjustable Shelf	2 pcs
Door	Double / Manual
External Size	1090x460x1650mm
Package Size	1210x580x1815 mm
Gross Weight	155 Kg

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com