

# WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET CBT11-4L



# WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET CBT11-4L

Weak Acid and Alkali Chemical Storage Cabinet are specifically designed for the safe storage of Weak Acids, Alkali, Hazardous Chemicals and corrosive materials for example Phosphoric acid, Carbonic acid, Acetic acid, Hydrocyanic acid and other chemical substances.

Used in Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants, Laboratory, Research. Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

## CBT11-4L WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET



### SPECIFICATIONS

Model	CBT11-4L
Capacity	4 gallon (15L)
Material	Double layer anti fire cold-rolled steel PP tray
Adjustable Shelf	1 pc
Door	Single / Manual
External Size	430x430x560 mm
Package Size	550x560x725 mm
Gross Weight	38 Kg

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@labstac.com](mailto:contact@labstac.com) | Website: [labstac.com](http://labstac.com)