

WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET CBT11-60L



WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET CBT11-60L

Weak Acid and Alkali Chemical Storage Cabinet are specifically designed for the safe storage of Weak Acids, Alkali, Hazardous Chemicals and corrosive materials for example Phosphoric acid, Carbonic acid, Acetic acid, Hydrocyanic acid and other chemical substances.

Used in Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants, Laboratory, Research. Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

CBT11-60L WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET



SPECIFICATIONS

Model	CBT11-60L
Capacity	60 gallon (227L)
Material	Double layer anti fire cold-rolled steel PP tray
Adjustable Shelf	2 pcs
Door	Double / Manual
External Size	860x860x1650 mm
Package Size	980x980x1815 mm
Gross Weight	195 Kg

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com