

# WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET CBT11-90L



# WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET CBT11-90L

Weak Acid and Alkali Chemical Storage Cabinet are specifically designed for the safe storage of Weak Acids, Alkali, Hazardous Chemicals and corrosive materials for example Phosphoric acid, Carbonic acid, Acetic acid, Hydrocyanic acid and other chemical substances.

Used in Laboratory, Research, Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants. Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

## CBT11-90L WEAK ACID AND ALKALI CHEMICAL STORAGE CABINET



### SPECIFICATIONS

Model	CBT11-90L
Capacity	90 gallon (340L)
Material	Double layer anti fire cold-rolled steel PP tray
Adjustable Shelf	2 pcs
Door	Double / Manual
External Size	1090x860x1650 mm
Package Size	1210x980x1815 mm
Gross Weight	243 Kg

**LABSTAC**

**Labstac LLC**

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: [contact@m.labstac.com](mailto:contact@m.labstac.com) | Website: [labstac.com](http://labstac.com)