

STRONG ACID AND ALKALI STORAGE CABINET CBT21-12L



STRONG ACID AND ALKALI STORAGE CABINET

CBT21-12L

Strong Acid and Alkali Storage Cabinet are specifically designed for the safe storage of Strong Acids, Alkali, Hazardous Chemicals and corrosive materials for example Hydrochloric acid, Sulfuric acid, nitric acid, KOH, NaOH and other chemical substances.

Used in Research, Pharmaceutical, Toxic Chemicals, Laboratory, Hazardous Contaminants, Potential Carcinogens. Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

CBT21-12L STRONG ACID AND ALKALI STORAGE CABINET

Padlock design, Adjustable shelves, Anti- corrosive materials for strong acid and alkali, Customizable



SPECIFICATIONS

| | |
|------------------|-----------------|
| Model | CBT21-12L |
| Capacity | 12 gallon (45L) |
| Material | pp |
| Adjustable Shelf | 1 pc |
| Door | Single / Manual |
| External Size | 590x460x890 mm |
| Package Size | 675x545x1025 mm |
| Gross Weight | 30 Kg |

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com