

LABSTAC

STRONG ACID AND ALKALI STORAGE CABINET CBT21-28L



labstac.com

STRONG ACID AND ALKALI STORAGE CABINET


CBT21-28L

Strong Acid and Alkali Storage Cabinet are specifically designed for the safe storage of Strong Acids, Alkali, Hazardous Chemicals and corrosive materials for example Hydrochloric acid, Sulfuric acid, nitric acid, KOH, NaOH and other chemical substances.

Used in Laboratory, Research, Pharmaceutical, Toxic Chemicals, Potential Carcinogens, Hazardous Contaminants. Also known as Acid Alkali Cabinet, Chemical Storage Cabinet.

CBT21-28L STRONG ACID AND ALKALI STORAGE CABINET

Padlock design, Adjustable shelves, Anti- corrosive materials for strong acid and alkali, Customizable



SPECIFICATIONS

Model	CBT21-28L
Capacity	28 gallon (105L)
Material	pp
Adjustable Shelf	2 pcs
Door	Double / Manual
External Size	910x600x900 mm
Package Size	995x685x1035 mm
Gross Weight	48 Kg

LABSTAC

Labstac LLC

82 Wendell Avenue, STE 100, Pittsfield, MA, 01201, USA

Email: contact@labstac.com | Website: labstac.com